

Product datasheet for **MR200957L3V**

Lyzl4 (NM_026915) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Lyzl4 (NM_026915) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Lyzl4
Synonyms:	1810009N24Rik; AI504271; Lyc4; Lyza
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_026915
ORF Size:	438 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR200957).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_026915.2
RefSeq Size:	1057 bp
RefSeq ORF:	438 bp
Locus ID:	69032
UniProt ID:	Q9D925
Cytogenetics:	9 F4
Gene Summary:	May be involved in fertilization (PubMed:21444326). Has no detectable bacteriolytic in vitro (PubMed:21444326). Has no lysozyme activity in vitro (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)