

Product datasheet for **MR200983L4V**

Magohb (NM_025564) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Magohb (NM_025564) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Magohb
Synonyms:	2010012C16Rik; Mago; Magoh; Magoh-rs1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_025564
ORF Size:	441 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR200983).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_025564.1 , NP_079840.1
RefSeq Size:	766 bp
RefSeq ORF:	447 bp
Locus ID:	66441
UniProt ID:	Q9CQL1
Cytogenetics:	6
Gene Summary:	Required for pre-mRNA splicing as component of the spliceosome. Plays a redundant role with MAGOH in the exon junction complex and in the nonsense-mediated decay (NMD) pathway.[UniProtKB/Swiss-Prot Function]



[View online »](#)