

Product datasheet for **MR200992A1V**

Mouse Ttr (NM_013697) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Mouse Ttr (NM_013697) AAV Particle
Tag: Myc-DDK
Symbol: Ttr
Synonyms: AA408768; AI787086; D17860; prea; prealbumin
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >MR200992 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGCTTCCTTCGACTCTTCCTCTTTGCCTCGCTGGACTGGTATTTGTGTCTGAAGCTGGCCCCGCGG
 GTGCTGGAGAATCCAAATGTCTCTGATGGTCAAAGTCCTGGATGCTGTCCGAGGCAGCCCTGCTGTAGA
 CGTGGCTGTAAAAGTGTTCAAAAGACCTCTGAGGGATCCTGGGAGCCCTTGCCTCTGGGAAGACCGCG
 GAGTCTGGAGAGCTGCACGGGCTCACCACAGATGAGAAGTTTGTAGAAGGAGGTACAGAGTAGAACTGG
 ACACCAATCGTACTGGAAGACACTTGGCATTTCCTCCGTTCCATGAATTCGCGGATGTGGTTTTCACAGC
 CAACGACTCTGGCCATCGCCACTACACCATCGCAGCCCTGCTCAGCCCATACTCCTACAGCACCACGGCT
 GTCGTCAGCAACCCCCAGAAT

ACGCGTACGCGGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR200992 protein sequence
 Red=Cloning site Green=Tags(s)

MASLRLFLLCLAGLVFVSEAGPAGAGESKCPLMVKVLDAVRGSPAVDVAVKVFKKTSEGSWEPFASGKTA
 ESGELHGLTTDEKFVEGVYRVELDTKSYWKTGLISPFHEFADVFTANDSGHRHYTIAALLSPYSYSTTA
 VVSNPQN

TRTRPLE**QKL**ISEEDLAANDILDYKDDDDKV

Species: Mouse
Serotype: AAV-2


[View online »](#)

ACCN:	NM_013697
ORF Size:	441 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_013697.1</u>
RefSeq Size:	1237 bp
RefSeq ORF:	444 bp
Locus ID:	22139
UniProt ID:	<u>P07309</u>
Cytogenetics:	18 11.47 cM
MW:	15.8 kDa

Product images:
