

Product datasheet for **MR201027L3V**

Ssbp1 (NM_212468) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Ssbp1 (NM_212468) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Ssbp1
Synonyms:	2810480P10Rik; G630031O20Rik; mtDBP; MtSSB
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_212468
ORF Size:	444 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR201027).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_212468.3 , NP_997633.1
RefSeq Size:	1489 bp
RefSeq ORF:	447 bp
Locus ID:	381760
UniProt ID:	Q9CYR0
Cytogenetics:	6 B1
Gene Summary:	This protein binds preferentially and cooperatively to ss-DNA. Probably involved in mitochondrial DNA replication. Associates with mitochondrial DNA (By similarity). [UniProtKB/Swiss-Prot Function]



[View online »](#)