

Product datasheet for MR201030

Crcp (NM_007761) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Crcp (NM_007761) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: Crcp
Synonyms: AL022669
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >MR201030 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGAAGTGAAGGATGCGAACGCAGCACTGCTCAGCAACTATGAGGTGTTCCAGTTACTAACGGATCTCA
 AAGAGCAGCGGAAGGAGAGTGGGAAGAACAAGCATAGCGCTGGGCAGCAGAACTGAACGCCATCACCTA
 TGAGACTAAAGTACATATCAAAAACACCATGTAGGAACCAGAGTCCAGCGATTGTCCAAGAATTCCTC
 ACAGCAATGAAAAGCCACAAGCTGACCAAAGCGGAGAAGCTGCAGCTGCTGAATCACAGGCCGATGACCG
 CGGTGGAGATCCAGCTGATGGTAGAAGAGAGTGAAGAACGGCTCACAGAGGAGCAGATCGAAGCCCTCCT
 GCACACTGTCACTAGCATTCTGCCTGCAGTCCAGAGGATGAACAGAGCAAAGCACTAGCAATGACGTG
 GCTATGGAAGAGGAGGAGCCAGCA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR201030 protein sequence
 Red=Cloning site Green=Tags(s)
 MEVKDANAALLSNYEVFQLLTDLKEQRKESGKNKHSAGQQLNLAITYETLKYISKTPCRNQSPAIVQEFLL
 TAMKSHKLTKAELQLLNHRPMTAVEIQLMVEESEERL TEEQIEALLHTVTSILPAGPEDEQSKSTSNDV
 AMEEEEPA

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI



[View online »](#)

Cloning Scheme:


ACCN: NM_007761

ORF Size: 447 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_007761.1](#), [NM_007761.2](#), [NP_031787.1](#)

RefSeq Size: 1474 bp

RefSeq ORF: 447 bp

Locus ID: 12909


UniProt ID: [O35427](#)

Cytogenetics: 5 G1.3

MW: 16.7 kDa

Gene Summary: DNA-dependent RNA polymerase catalyzes the transcription of DNA into RNA using the four ribonucleoside triphosphates as substrates. Specific peripheric component of RNA polymerase III which synthesizes small RNAs, such as 5S rRNA and tRNAs. Plays a key role in sensing and limiting infection by intracellular bacteria and DNA viruses. Acts as nuclear and cytosolic DNA sensor involved in innate immune response. Can sense non-self dsDNA that serves as template for transcription into dsRNA. The non-self RNA polymerase III transcripts induce type I interferon and NF- Kappa-B through the RIG-I pathway (By similarity). [UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR201030