

Product datasheet for **MR201036A1V**

Mouse (BC028925) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Mouse (BC028925) AAV Particle
Tag: Myc-DDK
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >MR201036 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGATGAGTCCTGCCAGTTCCTGTTTCTGCTAGTGCTCTGGATTCGGGAAACCAACGGTGATGTTGTGA
 TGACCCAGACTCCACTCACTTTGTCTGGTTACCATTGGACAACCAGCCTCCATCTCTTGAAGTCAAGTCA
 GAGCCTCTTAGATAGTGATGGAAAGACATATTTGAATTGGTTGTTACAGAGGCCAGGCCAGTCTCCAAAG
 CGCCTAATCTATCTGGTGTCTAACTGGACTCTGGAGTCCCTGACAGGTTCACTGGCAGTGGATCAGGGA
 CAGATTTACACTGAAAATCAGCAGAGTGGAGGCTGAGGATTTGGGAGTTTATTATTGCTGGCAAGGTAC
 ACATTTTCTGGACGTTCTGGTGGAGGCCAACAGCTGAAAATCAAACGGGCTGATGCTGCACCAACTGTAT
 CCATCTCCACCATCCAGTGAGCAGT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR201036 protein sequence
 Red=Cloning site Green=Tags(s)

MMSPAQFLFLLVLWIRETNQDVMTQTPLTSLVTIGQPASISCKSSQSLLDSDGKTYLNWLLQRPQGSPK
 RLIYLVSKLDSGVPDRFTGSGSGTDFTLKISRVEAEDLGVYYCQWGFHFGRSVEAPSWKSNGLMLHQLY
 PSSHPVSS

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Mouse
Serotype: AAV-2
ACCN: BC028925
ORF Size: 447 bp



[View online »](#)

Buffer: PBS with 0.001% Pluronic F68

Stability: AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.

RefSeq: [BC028925](#)

RefSeq Size: 1010 bp

RefSeq ORF: 449 bp

MW: 16.5 kDa