

## Product datasheet for MR201037

### Cend1 (NM\_021316) Mouse Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Cend1 (NM\_021316) Mouse Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** Cend1  
**Synonyms:** 1500001H12Rik; AI415214; BM88  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**ORF Nucleotide Sequence:** >MR201037 ORF sequence  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGAATCCCAGGAAAGTCAGCCAGCAGCCCCAAGCCAGATACCAAGGTGCCTCAGGCCACAGCTGAGG  
CCAAGGCCACTCCAGCTGCGGATGGGAAAGCTCCCTTGACCAAGCCTGTGAAGAAGGACACCCAAGCAGA  
GAAGCAGGAGCAAGCAGCAGCCCTGGGCCAGCCGCTACGAAGAAGACACCAGCCAAGGCAGATCCTGTA  
CTTCTCAACAACCACAGCAACCTGAAGCCAGCTCCACAGTGCCGGCAGCCCCAGCAGTCTGATGCAA  
CCTCAGAGCCCAAGGGTCTGGGGATGGGCAGAGGAGGACGAGTCCAACACTGGAGGCCGAGGTCCCTG  
GCCTTGTGAGAACTTGACCCCTTTGCTGGTAGCTGGGGGTGTGGCGTGGCCACCATAGCTTAATTCTT  
GGTGTGGCCTTTCTAGCCCGAAAAA

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >MR201037 protein sequence  
Red=Cloning site Green=Tags(s)

MESRGKSASSPKPDKVPQATAEAKATPAADGKAPLTKPVKDTQAEKQEQAAPGPAATKKTPAKADPV  
LLNNHNSNLKPAPTVPAPSSPDATSEPKGPGDGAEEDESNTGGRGPWCENLTPLLAVGGVAVATIALIL  
GVAFLARKK

**TR**TRPLEQKLISEEDLAANDILDYKDDDDKV

**Restriction Sites:** Sgfl-MluI



**Cloning Scheme:**


**ACCN:** NM\_021316

**ORF Size:** 450 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
  2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
  3. Close the tube and incubate for 10 minutes at room temperature.
  4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
  5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_021316.5](#)

**RefSeq Size:** 1657 bp


**RefSeq ORF:** 450 bp

**Locus ID:** 57754

**UniProt ID:** [Q9JKC6](#)

Cytogenetics: 7 F5  
 MW: 15 kDa  
 Gene Summary: Involved in neuronal differentiation.[UniProtKB/Swiss-Prot Function]

**Product images:**



Circular map for MR201037