

Product datasheet for MR201045

Cst6 (NM_028623) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Cst6 (NM_028623) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: Cst6
Synonyms: 1110017E11Rik; ichq; N28197
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >MR201045 representing NM_028623
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGAGCGTCCTCACTTCCCGCTGGCGATGGGCCTTGGCCTGCTCGATTCTGCCTCCTGACCCTGTAC
CTGACGCCCGCGCGAGCTGCGGAGCCGACGACCGGAGAACGCCAGAATTGTCACCCACCGACCCGCG
GGTGCAAAGGCCGCACAGGCAGCTGTGGCCAGCTACAATATGGGCAGCGACAGCCTCTACTACTCCGA
GACACCAAAGTCATCGATGCAAAATATCAGCTGGTGGCCTGATCAAGTACTACCTGACTTTGGACATAG
AAAGCACAGAGTCCGAAAGACCAGAGTCTCTGGAGAACACATGGACCTCACCACTGCCCTCTGGCTGC
AGGGGGGCAGCAGGAGAACTGCGTTGCAACTTTGAGCTCCTTGAAGTCCCTGGAAGAACAATACTCAG
CTTCTAAAGCATGACTGCGTGCAGGTG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR201045 representing NM_028623
Red=Cloning site Green=Tags(s)

MERPHFPLAMGLGLLAFCLLTLSPDARAE LRSRRTGERQNLSPDPRVQKAAQAAVASYNMGSDSLYYFR
DTKVIDAKYQLVAGIKYYLTLDIESTECKERTRVSGEHMDLTTCP LAAGGQQEKLRCNFELLEVPWKNTTQ
LLKHDCVQV

TRTRPLEQKLI SEEDLAANDILDYKDDDDKV


Chromatograms: https://cdn.origene.com/chromatograms/mm9005_c08.zip



[View online »](#)

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_028623

ORF Size: 447 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_028623.5](#), [NP_082899.1](#)

RefSeq Size: 3951 bp


RefSeq ORF: 450 bp

Locus ID: 73720

Cytogenetics: 19 4.3 cM

MW: 17.2 kDa

Product images:



Circular map for MR201045