

Product datasheet for **MR201064L4V**

Aif1l (NM_145144) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Aif1l (NM_145144) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Aif1l
Synonyms:	2810003C17Rik; AI043124; C87647; Iba2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_145144
ORF Size:	453 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR201064).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_145144.1 , NP_660126.1
RefSeq Size:	2967 bp
RefSeq ORF:	453 bp
Locus ID:	108897
UniProt ID:	Q9EQX4
Cytogenetics:	2 B
Gene Summary:	Actin-binding protein that promotes actin bundling. May neither bind calcium nor depend on calcium for function (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)