

Product datasheet for **MR201072L3V**

MyI6 (BC081470) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	MyI6 (BC081470) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	MyI6
Synonyms:	ESMLC, LC17A, LC17B, MLC1SM, MLC3SM, LC17-GI, LC17-NM, MLC3NM
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	BC081470
ORF Size:	453 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR201072).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	BC081470 , AAH81470
RefSeq Size:	646 bp
RefSeq ORF:	455 bp
Locus ID:	17904
Cytogenetics:	10 D3
Gene Summary:	Regulatory light chain of myosin. Does not bind calcium.[UniProtKB/Swiss-Prot Function]



[View online »](#)