

Product datasheet for **MR201075**

D16H22S680E (BC005445) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: D16H22S680E (BC005445) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: D16H22S680E
Synonyms: D16H22S680E; T10
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >MR201075 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGTGCATCATCTCTTTAAGTTTGATCCTCGTCCTGTTTCCAAAATGCATACAGGCTCATCTTGGCAG
CCAACAGAGATGAATTCTATAACAGGCCGTCGAAGTTGGCAGACTCTGGGGGAACAACAGTGAGATCCT
CAGTGGGCTCGACATGGAGGAAGGCAAGGCAGGAGGAACATGGCTGGGCATCAGCACACGTGGGAAGCTG
GGGGCACTCACCAACTACCTACAGCCCCGCAAGGCCAGACGCCGCTGGTAGAGGTGAACCTTGTGTCTC
ACTTTCTGACCAGTGACATGGACAGCCTGTCTACCTGAAGAAGGTCTCTACAGAGGGCCACCTGTATAA
CGGCTTTAACATCATAGCAGCTGACCTGAGCACGTCAAAGGAGATGTTGTTTGCTACTATGGAAACCGG
GGGGAGCCTGAGCCTATTGTCTTGACACCAGGC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR201075 protein sequence
Red=Cloning site Green=Tags(s)

MCIIFFKFDPRPVSKNAYRLILAAANRDEFYNRPSKLADFWGNSEILSGLDMEEGKAGGTWLGISTRGKL
GALTNYLQPRQEPDARGRGELVSHFLTSDMDSL SYLKKVSTEGHL YNGFNIIAADLSTSKGDVVCCYYGNR
GEPEPIVLTPG

TRTRPLEQKLI SEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN: BC005445

ORF Size: 453 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>BC005445</u> , <u>AAH05445</u>
RefSeq Size:	1607 bp
RefSeq ORF:	455 bp
Locus ID:	27883
Cytogenetics:	16 11.34 cM
MW:	16.7 kDa