

Product datasheet for **MR201086**

Fundc2 (NM_026126) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Fundc2 (NM_026126) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: Fundc2
Synonyms: 4833415N24Rik; DC44; HCBP6
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >MR201086 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGCCGCTAACAGTCAAGGAACTTTGATGGAAAATTTGAAGCACTGGATCTTGCAGAATTAATAAA
 AACAAACCATGGTGGCGTAAGCTGTTGGGCAGGAATCTGGGCCATCAGCTGAAAAATATAGTGTAGCAAC
 CCAGCTGGTAATTGGGGTGTACTGGATGGTGCCTGGTTTCGATTCCAGAAGTTGGAAAATTGGCT
 GCAACAGCTGTAGGAGGTGGATTCTTTCTCCTCAGCTGGCAAACCACTGGATACATCAAAGTTGACT
 GGCAAAGAGTAGAAAAAGACATGAAGAAAGCCAAAGAGCAGTTGAAAATTCGTAAGAATAAGCAAATACC
 AACTGAGGTCAAGAGCAAAGCTGAGGAGTTGTATCATTTGTAAGAAGAATGTTCTAGTACTGGAGGA
 TTTTTGGAGGCTTTCTGCTTGGCATGGCATCC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR201086 protein sequence
 Red=Cloning site Green=Tags(s)

MAANSQGNFDGKFEALDLAELTKKQPWWRKLFQESGPSAEKYSVATQLVIGGVGTWCTGFVQKVGKLA
 ATAVGGGFLLQLANHTGYIKVDWQRVEKDMKAKEQLKIRKNKQIPTEVKSAAEEVVSFVKKNVLVTGG
 FFGGFLMGAS

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: Sgfl-MluI



Cloning Scheme:


ACCN: NM_026126

ORF Size: 456 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_026126.4](#)

RefSeq Size: 3261 bp

RefSeq ORF: 456 bp


Locus ID: 67391

UniProt ID: [Q9D6K8](#)

Cytogenetics: X A7.3

MW: 16.6 kDa

Product images:



Circular map for MR201086