

Product datasheet for **MR201092A1V**

Mouse Ube2n (BC034898) AAV Particle

Product data:

Product Type: AAV Particles
 Product Name: Mouse Ube2n (BC034898) AAV Particle
 Tag: Myc-DDK
 Symbol: Ube2n
 Synonyms: 1500026J17Rik, UBC13
 Mammalian Cell Selection: None
 Vector: pAAV-AC-Myc-DDK (PS100089)
 ORF Nucleotide Sequence: >MR201092 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGC**

ATGGCCGGGCTGCCCCGAGGATCATCAAGGAAACCCAGCGTTTCTGGCAGAACCAAGTTCCTGGCATT
 AAGCAGAACCAAGATGAGAGCAACGCCGTTATTTTCATGTGGTCATTGCTGGCCCCAGGATTCCTCTT
 TGAGGGAGGGACTTTTAACTTGAAGTATTCCTCCAGAAGAATACCAATGGCAGCACCTAAAGTACGT
 TTCATGACCAAAATTTATCATCCTAATGTAGACAAGTTGGGAAGAATATGTTTAGATATTTGAAAGATA
 AGTGGTCCCAGCACTGCAGATCCGCACAGTTCTGCTATCAATCCAGGCTTTGTTAAGTGCTCCTAATCC
 AGATGATCCATTAGCAAATGATGTAGCCGAGCAATGGAAGACCAACGAAGCCCAAGCCATAGAAACAGCG
 AGAGCATGGACTAGGCTATATGCCATGAACAATATT

ACGCGTACGCGGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR201092 protein sequence
 Red=Cloning site Green=Tags(s)

MAGLPRIIKETQRLAEPVPGIKAEPDESNARYFHVVIAGPQDSPFEGGTFLKLEFLPEEYPMAPKVR
 FMTKIYHPNVDKLGRIKLDILKDKWSPALQIRTVLLSIQALLSAPNPDDPLANDVAEQWKTNEAQAIETA
 RAWTRLAMNNI

TRTRPLE**QKL**ISEEDLAANDILDYKDDDDKV

Species: Mouse
 Serotype: AAV-2



[View online »](#)

ACCN:	BC034898
ORF Size:	456 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>BC034898</u> , <u>AAH34898</u>
RefSeq Size:	743 bp
RefSeq ORF:	458 bp
Locus ID:	93765
Cytogenetics:	10 C2
MW:	17.1 kDa