

Product datasheet for **MR201096L3V**

Rps18 (NM_011296) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Rps18 (NM_011296) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Rps18
Synonyms:	H-2Ke3; H2-Ke3; Ke-3; ke3
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_011296
ORF Size:	459 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR201096).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_011296.1
RefSeq Size:	524 bp
RefSeq ORF:	459 bp
Locus ID:	20084
UniProt ID:	P62270
Cytogenetics:	17 17.98 cM
Gene Summary:	Located at the top of the head of the 40S subunit, it contacts several helices of the 18S rRNA. [UniProtKB/Swiss-Prot Function]


[View online »](#)