

## Product datasheet for MR201107

### Mgst3 (NM\_025569) Mouse Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Mgst3 (NM\_025569) Mouse Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** Mgst3  
**Synonyms:** 2010012L10Rik; 2010306B17Rik; 2700004G04Rik; AA516734; GST-III  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**ORF Nucleotide Sequence:** >MR201107 ORF sequence  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGATCGCC**

ATGGCTGTCTCTCTAAGGAGTATGGATTTGTGCTTCTCACTGGTGTGCCAGCTTTGTGATGGTGTCTCC  
 ACCTAGCCATCAACGTGGGCAAAGCCCGCAAGAAGTACAAGGTAGAGTACCCTGTCATGTACAGCACAGA  
 TCCTGAGAACGGGCATATGTTCAACTGCATTAGCGCGCCACCAGAACACGTTGGAGGTGTACCCTCCC  
 TTCTGTTTTTCTAACGGTGGGAGGTGTTTACCACCCGCGCATAGCTTCTGGCCTGGCCCTGGCCTGGA  
 TTATTGGGCGAGTCCTTTACGCATATGGCTACTACACAGGAGACCCTAGCAAGCGGTATCGAGGAGCCGT  
 GGGCTCTCTTGCCTCTTTGCCCTGATGGGCACCACCGTGTGCTCTGCTTCCAGCATCTCGGCTGGATC  
 AGACCAGGCTTAGGCTACGGGTCCAGATCCTGCCACCAC

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >MR201107 protein sequence  
 Red=Cloning site Green=Tags(s)  
 MAVLSKEYGFVLLTGAASFVMVLHLAINVGKARKKYKVEYPMYSTDPENGHMFNCIQRHQNTLEVYPP  
 FLFFLTVGGVYHPRIASGLGLAWIIGRVLYAYGYTGDPSKRYRGAVGSLALFALMGTTVCSAFQHLGWI  
 RPGLGYGSRSCHH

**TR**TRPLEQKLISEEDLAANDILDYKDDDDKV

**Restriction Sites:** Sgfl-MluI



[View online »](#)

**Cloning Scheme:**


**ACCN:** NM\_025569

**ORF Size:** 462 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_025569.2](#)

**RefSeq Size:** 635 bp


**RefSeq ORF:** 462 bp

**Locus ID:** 66447

**UniProt ID:** [Q9CPU4](#)

**Cytogenetics:** 1 H2.3  
**MW:** 17 kDa  
**Gene Summary:** Catalyzes oxydation of hydroxy-fatty acids. May participate to the lipid metabolism. [UniProtKB/Swiss-Prot Function]

**Product images:**



Circular map for MR201107