

Product datasheet for MR201125

Rnf122 (NM_175136) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Rnf122 (NM_175136) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: Rnf122
Synonyms: 1110063C11Rik
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >MR201125 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGCCACCATCCAGTGGTGAACGGGTGTTTCTGTGGCCTGGGACTGGTTAGTACCAACAAGTCCTGCT
CAATGCCACCGATCAGTTTCCAAGACCTCCCTCAATATCTACATGGTCATCTTCGGCACAGGCATCTT
CGTCTTCATGCTTAGCCTCATCTTCTGCTGCTACTTCATCAGCAAACCCGGAACCAGGCTCAGAGCGAG
CGATACGGCTACAAGGAGGTGGTCTTAAAGGTGATGCTAAGAAGTTACAGCTCTATGGACCTGTGCAG
TTTGCCTGGAAGACTTCAAGGGGAAGGATGAACTAGGCGTGCTTCCATGCCAGCATGCCTTTCACCGCAA
GTGCTGGTGAAGTGGCTAGAAGTACGTTGTGTCTGCCCATGTGCAACAAGCCATCGCGGCCCCACG
GAGACCTCGCAGAGCATCGGGATCCTGCTGGATGAACTGGTG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR201125 protein sequence
Red=Cloning site Green=Tags(s)
MHPFQWCNGCFGLGLVSTNKSCSMPPISFQDLPLNIYMVIFGTGIFVFMLSLIFCCYFISKLRNQAQSE
RYGYKEVVLKGDAAKLQLYGTCVLCLEDFKGDDELGVLPQHAHFHRKCLVKWLEVRCVCPMCNKPIAGPT
ETSQSIGILLDELV

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: Sgfl-MluI



Cloning Scheme:


ACCN: NM_175136

ORF Size: 465 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_175136.3](#)

RefSeq Size: 3556 bp


RefSeq ORF: 465 bp

Locus ID: 68867

UniProt ID: [Q8BP31](#)

Cytogenetics: 8 A3
 MW: 17.3 kDa
 Gene Summary: May induce necrosis and apoptosis. May play a role in cell viability (By similarity).
 [UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR201125