

## Product datasheet for **MR201140L3V**

### Ap1s3 (BC098237) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Ap1s3 (BC098237) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Ap1s3
Synonyms:	[s]3A, Jr2, 1190009B22
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	BC098237
ORF Size:	462 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR201140).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">BC098237.1</a>
RefSeq Size:	1298 bp
RefSeq ORF:	464 bp
Locus ID:	252903
Cytogenetics:	1 C4



[View online »](#)

**Gene Summary:**

Subunit of clathrin-associated adaptor protein complex 1 that plays a role in protein sorting in the late-Golgi/trans-Golgi network (TGN) and/or endosomes. The AP complexes mediate both the recruitment of clathrin to membranes and the recognition of sorting signals within the cytosolic tails of transmembrane cargo molecules (By similarity). Involved in TLR3 trafficking (By similarity).[UniProtKB/Swiss-Prot Function]