

## Product datasheet for MR201145L4V

## OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## Sod1 (NM\_011434) Mouse Tagged ORF Clone Lentiviral Particle

**Product data:** 

**Product Type:** Lentiviral Particles

**Product Name:** Sod1 (NM\_011434) Mouse Tagged ORF Clone Lentiviral Particle

Symbol: Sod

Synonyms: B430204E11Rik; Cu/Zn-SOD; CuZnSOD; Ipo-1; Ipo1; Sod-1; SODC

**Mammalian Cell** 

Selection:

Puromycin

**Vector:** pLenti-C-mGFP-P2A-Puro (PS100093)

Tag: mGFP

**ACCN:** NM\_011434

ORF Size: 465 bp

**ORF Nucleotide** 

The ORF insert of this clone is exactly the same as(MR201145).

Sequence:

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of

reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

**RefSeq:** <u>NM 011434.1</u>

 RefSeq Size:
 661 bp

 RefSeq ORF:
 465 bp

 Locus ID:
 20655

 UniProt ID:
 P08228

 Cytogenetics:
 16 51.56 cM

**Gene Summary:** Destroys radicals which are normally produced within the cells and which are toxic to

biological systems.[UniProtKB/Swiss-Prot Function]

