

Product datasheet for MR201153

Fundc1 (NM_028058) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Fundc1 (NM_028058) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: Fundc1
Synonyms: 1500005J14Rik; 1810033P05Rik
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >MR201153 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGCATCCCGGAACCCCTCCCAAGACTATGAAAGCGATGACGAATCATACGAAGTGTGGATTTAA
 CTGAATATGCAAGAAGACACCACTGGTGAATCGAGTATTTGGCCACAGTCCGGACCTATGGTAGAAAA
 ATACTCAGTAGCTACTCAGATTGTAATGGGTGGCGTGACTGGCTGGTGTGCAGGATTTTTATTCCAGAAG
 GTTGGAAAACCTGCTGCAACTGCAGTAGGTGGTGGTTTCCTTCTTTCAGGTTGCCAGTCACAGTGGCT
 ATGTACAGATCGACTGGAAGAGAGTTGAAAAGATGTAACAAAGCAAAGAGACAGATAAAGAAGCGAGC
 AAATAAAGCAGCACCTGAAATCAACAATAAATTGAAGAAGCAACAGACTTTATCAAGCAGAACATTGTG
 ATATCCAGCGGCTTCGTGGGAGGCTTTTTGCTAGGCCTGGCATCT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR201153 protein sequence
 Red=Cloning site Green=Tags(s)

MASRNPPPPQDYESDDESVEVLDL TEYARRHHWNRVFGHSSGPMVEKYSVATQIVMGGVTGWAGFLFQK
 VGKLAATAVGGGFLLLQVASHSGYVQIDWKRVEKDVNKAQRQIKKRANKAAPEINNIIEEATDFIKQNI
 ISSGFVGGFLLGLAS

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI



[View online »](#)

Cloning Scheme:


ACCN: NM_028058

ORF Size: 468 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_028058.4](#)

RefSeq Size: 1950 bp


RefSeq ORF: 468 bp

Locus ID: 72018

UniProt ID: [Q9DB70](#)

Cytogenetics: X A1.2
MW: 17.2 kDa
Gene Summary: Acts as an activator of hypoxia-induced mitophagy, an important mechanism for mitochondrial quality control.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR201153