

Product datasheet for **MR201189L4V**

Nudcd2 (NM_026023) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Nudcd2 (NM_026023) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Nudcd2
Synonyms:	2700047N05Rik; D11Ertd603e; RP23-382C18.4
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_026023
ORF Size:	474 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR201189).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_026023.4
RefSeq Size:	1599 bp
RefSeq ORF:	474 bp
Locus ID:	52653
UniProt ID:	Q9CQ48
Cytogenetics:	11 24.41 cM
Gene Summary:	May regulate the LIS1/dynein pathway by stabilizing LIS1 with Hsp90 chaperone. [UniProtKB/Swiss-Prot Function]



[View online »](#)