

Product datasheet for **MR201199L4V**

Bik (BC003732) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Bik (BC003732) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Bik
Synonyms:	Blk, Nbk
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	BC003732
ORF Size:	474 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR201199).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	BC003732 , AAH03732
RefSeq Size:	1379 bp
RefSeq ORF:	476 bp
Locus ID:	12124
Cytogenetics:	15 E1
Gene Summary:	Accelerates programmed cell death. Binding to the apoptosis repressors Bcl-X(L), BHRF1 or Bcl-2 suppresses this death-promoting activity.[UniProtKB/Swiss-Prot Function]



[View online »](#)