

Product datasheet for MR201205

Pgf (NM_008827) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids

Product Name: Pgf (NM_008827) Mouse Tagged ORF Clone

Tag: Myc-DDK

Symbol: Pgf

Synonyms: Al854365; PIGF; Plgf

Mammalian Cell Neomycin

Selection:

Vector:pCMV6-Entry (PS100001)E. coli Selection:Kanamycin (25 ug/mL)

ORF Nucleotide >MR201205 ORF sequence

Sequence: Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC

GCCGCGATCGCC

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT

ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR201205 protein sequence

Red=Cloning site Green=Tags(s)

MLVMKLFTCFLQVLAGLAVHSQGALSAGNNSTEVEVVPFNEVWGRSYCRPMEKLVYILDEYPDEVSHIFS PSCVLLSRCSGCCGDEGLHCVPIKTANITMQILKIPPNRDPHFYVEMTFSQDVLCECRPILETTKAERRK

TKGKRKRSRNSQTEEPHP

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: Sgfl-Mlul



OriGene Technologies, Inc. 9620 Medical Center Drive, Ste 200


CN: techsupport@origene.cn

Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com



Cloning Scheme:





^{*} The last codon before the Stop codon of the ORF

ACCN: NM_008827

ORF Size: 477 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of

reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube

containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method: 1. Centrifuge at 5,000xg for 5min.

2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.

3. Close the tube and incubate for 10 minutes at room temperature.

4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid

at the bottom.

5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of

shipping when stored at -20°C.

RefSeq: <u>NM 008827.3</u>

RefSeq Size: 1599 bp
RefSeq ORF: 477 bp
Locus ID: 18654
UniProt ID: P49764




Cytogenetics: 12 39.58 cM **MW:** 17.9 kDa

Gene Summary: Growth factor active in angiogenesis and endothelial cell growth, stimulating their

proliferation and migration. It binds to the receptor FLT1/VEGFR-1. Also promotes cell tumor

growth (By similarity).[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR201205