

Product datasheet for **MR201205L4V**

Pgf (NM_008827) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Pgf (NM_008827) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Pgf
Synonyms:	AI854365; PIGF; Plgf
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_008827
ORF Size:	477 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR201205).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_008827.2
RefSeq Size:	1599 bp
RefSeq ORF:	477 bp
Locus ID:	18654
UniProt ID:	P49764
Cytogenetics:	12 39.58 cM
Gene Summary:	Growth factor active in angiogenesis and endothelial cell growth, stimulating their proliferation and migration. It binds to the receptor FLT1/VEGFR-1. Also promotes cell tumor growth (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)