

Product datasheet for MR201208

4930505A04Rik (BC049719) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: 4930505A04Rik (BC049719) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: 4930505A04Rik
Synonyms: RP23-145N4.1
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >MR201208 representing BC049719
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGCAAAAGCCTTCTTGTGGAATAGTTCACCTACCAGTCTGGATGTATCAGGAGAACATGAGAACAATT
TTGTAGAATACATCTCTTTTATACATCAATATGATGCCAGGAGAACTCCCAATGAACCCATCAAAGGAAA
GAAACATGGAACCTTCGTGCAGAGAGAGATAAACTAGGGGCTATGCCAATAGTTCCTAAGGCACCGGAG
GTATTATTGAACACTCTAGAGTCAGTTTCATCAGAACAGCCCCAAAAACAGACAAAGGAAATTCATCA
GGAGACAAAGTGACCTCACCCAGGCTATGCCAACAGAATCCCAAGAGTTGCTGGAGACTTAGAACAGCT
TATCAGCACGAGAAGGCAGCCCAGCCCAGCAGACCCAGAAGCCTGGGTCTGGGAGATGACTGGAAGAGTG
TTGCTCTCAACTGGAGGGCAAGCTCTTCTCACCCAGGCACGGCCCTTAAACCAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA


Protein Sequence: >MR201208 representing BC049719
Red=Cloning site Green=Tags(s)
MQKPSGIVPLTSLDVSGEHENNFVEYISFIHQYDARRTPNEPIKGGKHGTFVQREIKLGAMPIVPKAPE
VLLNTLESGSSEQPKNRQRKFIRRQSDLTRPMPTEFPRVAGDLEQLISTRQSPADPEAWVWEMTGRV
LLSTGGQALLTRHGPLNH

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: Sgfl-MluI



Cloning Scheme:

Plasmid Map:


ACCN: BC049719

ORF Size: 474 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>BC049719</u>
RefSeq Size:	1081 bp
RefSeq ORF:	476 bp
Locus ID:	75087
Cytogenetics:	11 A4
MW:	39.6 kDa