

Product datasheet for MR201219L3

Hscb (NM_153571) Mouse Tagged Lenti ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Hscb (NM_153571) Mouse Tagged Lenti ORF Clone
Tag:	Myc-DDK
Symbol:	Hscb
Synonyms:	AI325508; AW049829; Hsc20
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR201219).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF.


ACCN:	NM_153571
ORF Size:	477 bp



[View online »](#)

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	NM_153571.1 , NP_705799.1
RefSeq Size:	825 bp
RefSeq ORF:	705 bp
Locus ID:	100900
UniProt ID:	Q8K3A0
Cytogenetics:	5 F
Gene Summary:	Acts as a co-chaperone in iron-sulfur cluster assembly in both mitochondria and the cytoplasm. Required for incorporation of iron-sulfur clusters into SDHB, the iron-sulfur protein subunit of succinate dehydrogenase that is involved in complex II of the mitochondrial electron transport chain. Recruited to SDHB by interaction with SDHAF1 which first binds SDHB and then recruits the iron-sulfur transfer complex formed by HSC20, HSPA9 and ISCU through direct binding to HSC20. Also mediates complex formation between components of the cytosolic iron-sulfur biogenesis pathway and the CIA targeting complex composed of CIAO1, DIPK1B/FAM69B and MMS19 by binding directly to the scaffold protein ISCU and to CIAO1. This facilitates iron-sulfur cluster insertion into a number of cytoplasmic and nuclear proteins including POLD1, ELP3, DPYD and PPAT.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR201219L3