

## Product datasheet for **MR201227L3V**

### IL1rn (NM\_031167) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	IL1rn (NM_031167) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	IL1rn
Synonyms:	F630041P17Rik; IL-1ra
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_031167
ORF Size:	480 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR201227).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_031167.5</a> , <a href="#">NP_112444.1</a>
RefSeq Size:	2464 bp
RefSeq ORF:	480 bp
Locus ID:	16181
UniProt ID:	<a href="#">P25085</a>
Cytogenetics:	2 16.36 cM
Gene Summary:	Inhibits the activity of interleukin-1 by binding to receptor IL1R1 and preventing its association with the coreceptor IL1RAP for signaling. Has no interleukin-1 like activity. [UniProtKB/Swiss-Prot Function]



[View online »](#)