

Product datasheet for **MR201255L4V**

Taf12 (NM_025579) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Taf12 (NM_025579) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Taf12
Synonyms:	20kDa; 2810422D08Rik; AW557038; Taf2]
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_025579
ORF Size:	486 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR201255).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_025579.1
RefSeq Size:	1353 bp
RefSeq ORF:	486 bp
Locus ID:	66464
UniProt ID:	Q8VE65
Cytogenetics:	4 D2.3



[View online »](#)

Gene Summary:

TAFs are components of the transcription factor IID (TFIID) complex, PCAF histone acetylase complex and TBP-free TAFII complex (TFTC). TAFs components-TIIFD are essential for mediating regulation of RNA polymerase transcription (By similarity).[UniProtKB/Swiss-Prot Function]