

Product datasheet for **MR201272**

6330527O06Rik (BC004791) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	Myc-DDK
Symbol:	6330527O06Rik
Synonyms:	3110035N03Rik; 6330527O06Rik; BAD-LAMP; LAMP-5
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)

ORF Nucleotide Sequence: >MR201272 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGCTCTTTGTAAGGAAAGTCACAACACTTCCAAAGGACCAGAGGCGACTTGAATCTGAACAAGGTGC
ATTTTGTCTATGATTCTCGGAGAAGACTCACTTTAAAGCCCTGTGAAGGTTAATAAGTACATCGCCAG
CTCTCATCACCTCTCTGCCTTGGTCACCCAGCTGGGATGTCCTATGAGTGTGAGGCTCAGCAGACCATT
TCTCTGGCCTCCAGTGACCCCTCAGAAGACTGTCAACATGATCCTGTCCGAGTGCACATCCAGCCCTTTG
ATATCATCTCTGACTTTTGTCTTCAGTGAAGAACATAAATGTCCAGTGGATGAGCAAGAGCAGTTGGAAGA
GACCCCTGCCCTGATCCTGGGCCTCATCTTGGGCCTTGTATTGTCATAACACTTGTGATCTACCACATC
CACCATAAAATGACTGCCAATCAAGTGCAGATCCCCAGGGACCGTTCCAGTACAAGCACATGGGC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR201272 protein sequence
Red=Cloning site Green=Tags(s)

MLFVKESHNTSKGPEATWNLNKVHFVYDSSEKTHFKAPVKVNKYIASSHLSALVTPAGMSYECQAQQT
SLASSDPQKVTMILSAVHIQPFDIISDFVFSSEHKCPVDEQEQLLEETLPLILGLILGLVIVITLVIYHI
HHKMTANQVQIPDRSQYKHM

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI



Cloning Scheme:



ACCN: BC004791

ORF Size: 486 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:


1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [BC004791](#), [AAH04791](#)

RefSeq Size:	1237 bp
RefSeq ORF:	488 bp
Locus ID:	76161
Cytogenetics:	2 F3
MW:	18.3 kDa
Gene Summary:	Plays a role in short-term synaptic plasticity in a subset of GABAergic neurons in the brain. [UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR201272