

Product datasheet for **MR201276L3V**

Cav2 (NM_016900) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Cav2 (NM_016900) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Cav2
Synonyms:	AI447843
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_016900
ORF Size:	486 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR201276).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_016900.4 , NP_058596.1
RefSeq Size:	2748 bp
RefSeq ORF:	489 bp
Locus ID:	12390
UniProt ID:	Q9WVC3
Cytogenetics:	6 A2



[View online »](#)

Gene Summary:

This gene belongs to the caveolin family whose members encode the major protein components of caveolae, which are invaginations of plasma membrane. This gene is located adjacent to caveolin-1 and the proteins coexpressed by the two genes localize together in caveolae, where they form hetero-oligomers. The encoded protein may be involved in diverse cellular functions including proliferation, differentiation, endocytosis and trafficking. Alternative splicing of this gene results in transcript variants encoding different isoforms. [provided by RefSeq, Apr 2013]