

Product datasheet for **MR201281L4V**

Ly86 (NM_010745) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Ly86 (NM_010745) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Ly86
Synonyms:	ly86_tv2; MD-1; MD1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_010745
ORF Size:	489 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR201281).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_010745.1 , NP_034875.1
RefSeq Size:	997 bp
RefSeq ORF:	489 bp
Locus ID:	17084
UniProt ID:	O88188
Cytogenetics:	13 A3.3
Gene Summary:	May cooperate with CD180 and TLR4 to mediate the innate immune response to bacterial lipopolysaccharide (LPS) and cytokine production. Important for efficient CD180 cell surface expression.[UniProtKB/Swiss-Prot Function]



[View online »](#)