

Product datasheet for MR201286

Lysmd4 (BC057459) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Lysmd4 (BC057459) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: Lysmd4
Synonyms: 4930506D23Rik
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >MR201286 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGCCCTTGGGCTTCATCCTCTGAGACCAGGGTGACCCTGGTAGACCTGCCTGAAGATGAAGATGCTG
 GTGGAGCAACCACTCAGGGCAACCAGCTGACAGACTTCTCAAGGGCATTGATGAGAACATTGAGCGGGC
 TGTGCATTCTGACGTCTCCACGGTGACAGCTGCTGTGTGGAGGCCCCAGATCAGCTGCTGCTCCCCATA
 ACCCAGAAGCCAGTAGCCGATGGTGCCGACTGTGGGATTCAGTGGTGAATGCCGTCTTCCTTATGCTTC
 TCATTGGGATCGTATTGCCAGTGTCTATCTCGTCTACTTCAAAATACAGGCCACCGGGGAGCCCTCAA
 TGGCTTGAACGCAACTGTTGTCCCAATGGCTCGATGACACTGAGCCCAGTTCAGGGCAGGCTCCAGG
 TTGGCAATCCCAGTGCCTACCCTCCCCGCTTCAGACAGCCAGGTACGCCCTACCACCCAGGCAGGGGCC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR201286 protein sequence
 Red=Cloning site Green=Tags(s)

MPLGASSETRVTLVDLPEDEDAGGATTQGNQLTDFFKGIDENIERAVHSDVDFHGDSCCVEAPDQLLLPI
 TQKPVADGADCGIQWNAVFLMLLIGIVLPVFLVYFKIQATGEPNSGLNATVVPNGSMTLSPVPGQAPR
 LAIPVPTLPASDSQVSPPTQAGA

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI



Cloning Scheme:


ACCN: BC057459

ORF Size: 489 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [BC057459](#), [AAH57459](#)

RefSeq Size: 2954 bp


RefSeq ORF: 491 bp

Locus ID: 75099

Cytogenetics: 7 C

MW: 17.1 kDa

Product images:



Circular map for MR201286