

Product datasheet for MR201299

Emp3 (NM_001146346) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Emp3 (NM_001146346) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: Emp3
Synonyms: H-4; H4; HNMP-1; MI-35; Ymp
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >MR201299 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGTCACTCCTCCTGTTGGTGGTCTCTGCCCTCCACATCCTCATTCTTGCTTGCTTTTTGTGGCCACTT
TGGACAAGTCCTGGTGGACTCTCCAGACAAAGAGTCCCTGAACCTGTGGTATGACTGCACGTGGAACAC
CACCACCAAACATGGGCTGCAGTAACGTCAGTGAGAATGGCTGGCTGAAGGCAGTGCAGGCGCTCATG
GTGCTGTCTCTCATCCTCTGCTGCCTTCCTCCTCCTCATGTTCCAACCTACACCATGCGGCGCG
GAGGGCTCTTCTACGCTACCGGCTCTGCCAGCTTGCACCAGTGCAGCTGTGTTCTCCGGGCACTCAT
CTATGCCATCCACACCGAGGAGATCCTGGCCAAGCACCCGAGTGGGGGCACTTCGGTACTGCTTCGCC
CTGGCCTGGGTGGCTTTTCCACTCGCTCTGGTCAGCGGCATTGTCTACATCCACCTGCGGAAACGTGAA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR201299 protein sequence
Red=Cloning site Green=Tags(s)

MSLLLLVVSALHILILVLLFVATLTKSWWTLDPKESLNLWYDCTWNTTQTWACSNVSENGWLKAVQALM
VLSLILCCLSFILFMFQLYTMRRGGLFYATGLCQLCTSAAVFSGALIYAIHTEEILAKHPSGGSFGYCF
LAWVAFPLALVSGIVYIHLRKRE

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI



[View online >](#)

Cloning Scheme:


ACCN: NM_001146346

ORF Size: 492 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001146346.1](#), [NP_001139818.1](#)

RefSeq Size: 787 bp


RefSeq ORF: 492 bp

Locus ID: 13732

UniProt ID: [O35912](#)

Cytogenetics: 7 29.6 cM
MW: 18.3 kDa
Gene Summary: Probably involved in cell proliferation and cell-cell interactions.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR201299