

## Product datasheet for **MR201304L4V**

### BC003885 (BC003885) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	BC003885 (BC003885) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	BC003885
Synonyms:	MGC6735
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	BC003885
ORF Size:	489 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR201304).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">BC003885</a> , <a href="#">AAH03885</a>
RefSeq Size:	1568 bp
RefSeq ORF:	491 bp
Locus ID:	225215
Cytogenetics:	9 D
Gene Summary:	Involved in the biogenesis of the 60S ribosomal subunit. Ensures the docking of GTPBP4/NOG1 to pre-60S particles (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)