

Product datasheet for **MR201331**

Tceal1 (NM_146236) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Tceal1 (NM_146236) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: Tceal1
Synonyms: 0610011M09Rik; BC011290; P21; pp21; SIIR
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >MR201331 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGAAAAACACGCGCAGTGAAAAAGAAGAGCAACCAGAGAGCACTCTGAAGATCGATGAAGAGCAGCCTG
CAGTGGAGCAGTCTCCTGAAAAATCAGTGCTCTGAGGAGGATCAATCCTCGGAGGATCTGTCTCCGAGGA
GCAGTCTCCGAGGAGGAGTTCTTTCTGAAGAGCTCTGCCGAGCTGCTCCAGAGATGCTGCTGTCT
GAAGATCGTCTCCGAGGAGTGCCTCACAGAAGAACCAGTTCGAAGACCGTATCCCATGGAGCAAC
CTCCCTGCGGGTGGGCAACACAAGCTAGAGGAGGGAAGCTTCAAAGAAAGACTGGCCCGTATCCGCC
TCAATTCATAGGAGACATTCATGGCAGAAATTTAAGTAATGAAGAAATGATACAGGCAGCGGATGAGTTG
GAAGAGATGAAGAGAGTACGAAACAACTGATGATAATGCATTGGAAGGCAAACCGAGCCGCCCTTACC
CTATT

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR201331 protein sequence
Red=Cloning site Green=Tags(s)

MENTRSENEEQPESTLKIIEEQPAVEQSPENQCSEEDQSSSEDLSSSEEQSSEEEFFPEELLPELLPEMLLS
EDRPPQECLSQKNQFEDRIPMEQPCCGVGKHKLEEGSFKERLARIRPQFIGDIHGRNLSNEEMIQADEL
EEMKRVNRKLMIMHWKAKRSRPYPI

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: Sgfl-MluI



[View online »](#)

Cloning Scheme:


ACCN: NM_146236

ORF Size: 498 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_146236.2](#)

RefSeq Size: 1124 bp

RefSeq ORF: 498 bp

Locus ID: 237052


UniProt ID: [Q921P9](#)

Cytogenetics: X F1

MW: 19.3 kDa

Gene Summary: May be involved in transcriptional regulation. Modulates various viral and cellular promoters in a promoter context-dependent manner. Does not bind DNA directly (By similarity). [UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR201331