

## Product datasheet for **MR201331L4V**

### Tceal1 (NM\_146236) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Tceal1 (NM_146236) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Tceal1
Synonyms:	0610011M09Rik; BC011290; P21; pp21; SIIR
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_146236
ORF Size:	498 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR201331).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_146236.1</a>
RefSeq Size:	1124 bp
RefSeq ORF:	498 bp
Locus ID:	237052
UniProt ID:	<a href="#">Q921P9</a>
Cytogenetics:	X F1
Gene Summary:	May be involved in transcriptional regulation. Modulates various viral and cellular promoters in a promoter context-dependent manner. Does not bind DNA directly (By similarity). [UniProtKB/Swiss-Prot Function]



[View online »](#)