

Product datasheet for MR201335L3

Pmch (NM_029971) Mouse Tagged Lenti ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Pmch (NM_029971) Mouse Tagged Lenti ORF Clone
Tag:	Myc-DDK
Symbol:	Pmch
Synonyms:	A230109K23Rik
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR201335).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF.


ACCN:	NM_029971
ORF Size:	498 bp



[View online »](#)

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_029971.2 , NP_084247.1
RefSeq Size:	748 bp
RefSeq ORF:	498 bp
Locus ID:	110312
UniProt ID:	P56942
Cytogenetics:	10 43.7 cM
Gene Summary:	This gene encodes a preproprotein that is cleaved to produce three mature peptides: melanin concentrating hormone, neuropeptide-glutamic acid-isoleucine (NEI), and neuropeptide-glycine-glutamic acid (NGE). The proprotein is processed differently depending on the tissue where it is expressed. Melanin concentrating hormone is involved in the regulation of food intake, energy homeostasis, and sleep-wake behavior. Disruption of this gene is associated with resistance to diet-induced obesity. [provided by RefSeq, May 2013]

Product images:



Circular map for MR201335L3