

Product datasheet for **MR201341L4V**

Hint3 (NM_025798) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Hint3 (NM_025798) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Hint3
Synonyms:	0610010I17Rik; AV009015; HINT-3; HINT-4; HINT4
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_025798
ORF Size:	498 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR201341).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_025798.2
RefSeq Size:	696 bp
RefSeq ORF:	498 bp
Locus ID:	66847
UniProt ID:	Q9CPS6
Cytogenetics:	10 A4
Gene Summary:	Hydrolyzes phosphoramidate and acyl-adenylate substrates.[UniProtKB/Swiss-Prot Function]



[View online »](#)