

## Product datasheet for **MR201356L3V**

### Mtfp1 (NM\_026443) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Mtfp1 (NM_026443) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Mtfp1
Synonyms:	1700020C11Rik; 2610507A21Rik; AV005788; Mtp18
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_026443
ORF Size:	498 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR201356).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_026443.3</a>
RefSeq Size:	1397 bp
RefSeq ORF:	501 bp
Locus ID:	67900
UniProt ID:	<a href="#">Q9CRB8</a>
Cytogenetics:	11 A1
Gene Summary:	Involved in the mitochondrial division probably by regulating membrane fission. Loss-of-function leads to apoptosis (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)