

Product datasheet for **MR201357L4V**

Dgcr6 (BC048503) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Dgcr6 (BC048503) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Dgcr6
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	BC048503
ORF Size:	498 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR201357).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	BC048503.1
RefSeq Size:	948 bp
RefSeq ORF:	500 bp
Locus ID:	13353
Cytogenetics:	16 11.19 cM



[View online »](#)

Gene Summary:

This gene encodes a protein that is similar to the gonadal protein in *Drosophila* (fruit fly). The encoded protein is thought to play a role in migration of neural crest cells during development. Deletions in the human gene are associated with DiGeorge syndrome (or velocardiofacial syndrome) which has many clinical features including cardiac abnormalities, cleft palate, atypical facial features, hypocalcemia, hypoparathyroidism and defective development or congenital absence of the thymus. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jan 2014]