

Product datasheet for **MR201361**

Cxcl15 (NM_011339) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Cxcl15 (NM_011339) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Cxcl15
Synonyms:	Il8; lungkine; Scyb15; weche
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>MR201361 representing NM_011339 Red=Cloning site Blue=ORF Green=Tags(s)

TTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCTGCTCAAGGCTGGTCCATGCTCCTGCTGGCTGTCCTTAACCTAGGCATCTTCGTCCTGCCCTGTG
ACACTCAAGAGCTACGATGTCTGTGATTACAGGAACACTCTGAATTCATTCCTCTCAAACCTATTA
TATAATGGTGATATTCGAGACCATTTACTGCAACAGAAAGGAAGTGATAGCAGTCCCAAAAAATGGGAGT
ATGATTTGTTGGATCCTGATGCTCCATGGGTGAAGGCTACTGTTGGCCCAATTACTAACAGGTTCTAC
CTGAGGACCTCAAACAAAGGAATTTCCACCGCAATGAAGCTTCTGTATAGTGTGAGCATGAAAAGCC
TCTATATCTTTTCAATTTGGGAGACCTGAGAACAGAGAATATTTCCCTTTCCAATTCGGGAGACCTCTAGA
CACTTTGCTGATTTAGCTCACAACAGTGATAGGAATTTTCTACGGGACTCCAGTGAAGTCAGCTTGACAG
GCAGTGATGCC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Chromatograms: https://cdn.origene.com/chromatograms/mm9050_e03.zip

Restriction Sites: SgfI-MluI



[View online »](#)

Cloning Scheme:



ACCN: NM 011339

ORF Size: 501 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:


1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: NM_011339.2, NP_035469.1

RefSeq Size: 2127 bp

RefSeq ORF:	504 bp
Locus ID:	20309
UniProt ID:	Q9WVL7
Cytogenetics:	5 44.78 cM
MW:	19.5 kDa
Gene Summary:	Chemotactic for neutrophils. Involved in lung-specific neutrophil trafficking during normal and inflammatory conditions.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR201361