

Product datasheet for **MR201365L4V**

Cetn3 (NM_007684) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Cetn3 (NM_007684) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Cetn3
Synonyms:	MmCEN3
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_007684
ORF Size:	504 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR201365).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_007684.3 , NP_031710.1
RefSeq Size:	1013 bp
RefSeq ORF:	504 bp
Locus ID:	12626
UniProt ID:	O35648
Cytogenetics:	13 C3
Gene Summary:	Plays a fundamental role in microtubule-organizing center structure and function. [UniProtKB/Swiss-Prot Function]



[View online »](#)