

Product datasheet for **MR201375L3V**

Cggbp1 (NM_178647) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Cggbp1 (NM_178647) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Cggbp1
Synonyms:	AA960172; AL023996; AW045841
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_178647
ORF Size:	504 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR201375).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_178647.2
RefSeq Size:	4382 bp
RefSeq ORF:	504 bp
Locus ID:	106143
UniProt ID:	Q8BHG9
Cytogenetics:	16 C1.3
Gene Summary:	Binds to nonmethylated 5'-d(CGG)(n)-3' trinucleotide repeats in the FMR1 promoter. May play a role in regulating FMR1 promoter (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)