

## Product datasheet for MR201397

### Cd302 (BC005501) Mouse Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Cd302 (BC005501) Mouse Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** Cd302  
**Synonyms:** 1110055L24Rik; A1159627  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**ORF Nucleotide Sequence:** >MR201397 ORF sequence  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCCCGCATCGCC

ATGGTAAGCATACACAATGAAGAGGAAAACGCGTTTATACTGGACACTTTGCAAAGCGATGGAAGGGTC  
 CAGATGATCTCCTGCTAGGCATGTTCTATGACACTGATGATGCAACTTCAAGTGGTATGATCATTCAA  
 TATGACATTCGACAAGTGGCAGATCAAGATGGTGAGGACCTAGTTGATACCTGTGGTTTTCTGTACACC  
 AAGACAGGTGAATGGAGAAAAGGGGATTGTGAAATCTCTCTGTGGAGGGAACACTTTGCAAAGCAGCAA  
 TCCCATATGACAAGAAGTATTTATCAGATAACCACATTTAATATCGACTCTGGTGATCGCTAGCACAGT  
 AACTCTGGCAGTTTTGGGAGCGATCATTGGTTCTCTATAGAAGAAACGCGCGCTCTGGCTTACCTCT  
 TTTTCACCTGCACCACTGTCACCTTACAGTGATGGCTGTGCCCTGGTAGTTGCAGAAGAAGATGAATATG  
 CTGTTCAGCTGGAC

AGCGGACCGACGCGTACGCGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC  
 TGGATTACAAGGATGACGACGATAAGGTTTAA


**Protein Sequence:** >MR201397 protein sequence  
 Red=Cloning site Green=Tags(s)

MVSIHNEEENAFILDTLQKRWGPDDLILLGMFYDTHDQDFKWDHNSMTFDKWADQDGEDLVDTCGFLYT  
 KTGEWRKGDCEISSVEGTLCKAAIPYDKKYLSDNHILISTLVIASVTVLAVLGAIWFLYRRNARSGFTS  
 FSPAPLSPYSDGCALVVAEEDYAVQLD

SGPTRTRPLEQKLI SEEDLAANDILDYKDDDDKV

**Restriction Sites:** Sgfl-RsrII



**Cloning Scheme:**


**ACCN:** BC005501

**ORF Size:** 504 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [BC005501](#), [AAH05501](#)

**RefSeq Size:** 798 bp

**RefSeq ORF:** 506 bp


**Locus ID:** 66205

**Cytogenetics:** 2 C1.1

**MW:** 19 kDa

**Gene Summary:** Potential multifunctional C-type lectin receptor that may play roles in endocytosis and phagocytosis as well as in cell adhesion and migration.[UniProtKB/Swiss-Prot Function]

**Product images:**



Circular map for MR201397