

Product datasheet for MR201398

Dazap2 (NM_011873) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Dazap2 (NM_011873) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: Dazap2
Synonyms: AI314727; Brbp; Gcap28; gt6-12; mKIAA0058; Prtb
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >MR201398 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGAACAGCAAAGGTCAATATCCAACACAGCCTACCTACCCTGTGCAACCTCCTGGGAATCCAGTATATC
 CTCAGACCTTACATCTTCCTCAGGCTCCACCCTATACTGATGCTCCACCTGCATATTCAGAGCTCTATCG
 TCCAAGCTTTGTTTCATCCAGGGGCTGCCACAGTCCCCACCATGTCAGCTGCATTTCTGGCGCTCATTG
 TATCTGCCTATGGCCAGTCTGTGGCTGTTGGCCTCTAGGCTCCACAATCCCCATGGCTTATTATCCAG
 TTGGTCCCATCTATCCACCTGGCTCAGCAGTGTGGTGAAGGAGGGTATGATGCAGGTGCCAGATTTGG
 AGCTGGTGCCACTGCTGGCAACATTCTCCTCCACCTCCTGGATGCCCTCCCAACGCTGCTCAGCTCGCA
 GTCATGCAGGGAGCCAACGTCCTTGAACCTCAGCGCAAGGAACTTCTTCATGGGTGGCTCAGATGGTG
 GCTACACCATTGG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR201398 protein sequence
 Red=Cloning site Green=Tags(s)
 MNSKGQYPTQPTYVPVQPPGNPVYPQTLHLPQAPPYTDAPPAYSELYRPSFVHPGAATVPTMSAAFPGASL
 YLPMAQSVAVGPLGSTIPMAYYPVGPIYPPGSAVLVEGGYDAGARFGAGATAGNIPIPPPPGCPPNAAQLA
 VMQGANVLVTQRKGNFFMGGSDGGYTIW

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: Sgfl-MluI



Cloning Scheme:


ACCN: NM_011873

ORF Size: 507 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_011873.2](#), [NP_036003.2](#)

RefSeq Size: 1860 bp

RefSeq ORF: 507 bp


Locus ID: 23994

UniProt ID: [Q9DCP9](#)

Cytogenetics: 15 56.36 cM

MW: 17.3 kDa

Product images:



Circular map for MR201398