

Product datasheet for MR201412L2V

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Cdkn2a (NM 009877) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type: Lentiviral Particles

Product Name: Cdkn2a (NM_009877) Mouse Tagged ORF Clone Lentiviral Particle

Symbol:

Arf; ARF-INK4a; INK4a-ARF; Ink4a/Arf; MTS1; p16; p16(INK4a); p16INK4a; p19 p19ARF; Pctr1 Synonyms:

Mammalian Cell

Selection:

None

pLenti-C-mGFP (PS100071) Vector:

mGFP Tag:

ACCN: NM 009877

ORF Size: 507 bp

OTI Disclaimer:

ORF Nucleotide

Sequence:

The ORF insert of this clone is exactly the same as(MR201412).

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This

clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

RefSeq: NM 009877.2, NP 034007.1

RefSeq Size: 929 bp RefSeq ORF: 510 bp 12578 Locus ID: **UniProt ID:** Q64364

Cytogenetics: 4 42.15 cM

Gene Summary: Acts as a negative regulator of the proliferation of normal cells by interacting strongly with

CDK4 and CDK6. This inhibits their ability to interact with cyclins D and to phosphorylate the

retinoblastoma protein.[UniProtKB/Swiss-Prot Function]

