

Product datasheet for **MR201430L3V**

Btg1 (BC006834) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Btg1 (BC006834) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Btg1
Synonyms:	AI426953; AW546738
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	BC006834
ORF Size:	513 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR201430).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	BC006834 , AAH06834
RefSeq Size:	1830 bp
RefSeq ORF:	515 bp
Locus ID:	12226
Cytogenetics:	10 C3
Gene Summary:	Anti-proliferative protein.[UniProtKB/Swiss-Prot Function]



[View online »](#)