

Product datasheet for MR201438

Nsg2 (NM_008741) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Nsg2 (NM_008741) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: Nsg2
Synonyms: 8.5; AA989750; P19; R75287
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >MR201438 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCCGCATCGCC

ATGGTGAAGCTGAACAGCAACCCTGGAGAGAAGGGGGCCAAGCCGCCTTCGGTGGAGGATGGTTTCCAGA
 CCGTCCCTCTCATCAGCCCTTGGAGGTTAATCATTTGCAGCTGGCTGCTCCGAGAAGGTGATCGTGAA
 GACAAGAACAGAATACCAGCCGGAACAGAGGAACAAAGGGAAGTCCGGGTGCCAAGATCGCTGAATTT
 ACGGTACCAATTCTGTGAGCCTGGCCCTTGCTTCCTTGCTTGATTGTGTTCTGGTGGTTTACAAAG
 CCTTTACCTATGACCACAGCTGTCCTGAAGGATTTGTCTACAAGCACAAGCGCTGTATCCCGGCCTCGCT
 GGATGCTTACTACTCGTCCCAGGACCCAGCTCCAGGAGCCGCTTCTACACAGTCATCAGCCACTACAGT
 GTGGCCAAGCAGAGCACTGCCCGGGCCATCGGGCCATGGCTGTCAGCAGCTGCTGTCATCCATGAGCCCA
 AGCCACCCAAGACGCAGGGCCAT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR201438 protein sequence
 Red=Cloning site Green=Tags(s)
 MYKLNSNPGEKGAKPPSVEDGFQTVPLITPLEVNLQLAAPEKVIKTRTEYQPEQRNKGKFRVPKIAEF
 TVTILVSLALAFACIVFLVYKAFETYDHSCEPGEFVYKHKRCIPASLDAYSSQDPSSRSRYFTVISHYS
 VAKQSTARAI GPWLSAAAVIHEPKPKTQGH

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI


[View online »](#)

Cloning Scheme:



ACCN: NM_008741

ORF Size: 513 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_008741.4](#)

RefSeq Size: 2308 bp

RefSeq ORF: 516 bp


Locus ID: 18197

UniProt ID: [P47759](#)

Cytogenetics: 11 A4

MW: 19 kDa

Product images:



Circular map for MR201438