

Product datasheet for **MR201446**

Spata3 (NM_027029) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Spata3 (NM_027029) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: Spata3
Synonyms: 1700011N12Rik; 1700029H01Rik; 4930424D10Rik; Mtsarg1; Tsarg1
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >MR201446 representing NM_027029
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGAAGAAGGTCAAAAAGAAAAGTCTGATTCCAGGCGCCGAGGAATTCATTTACCCGAGACTAGCT
CCGATTCCTCACAGCAGCCAGCTCAGAAACCCGCGTCATGTCCAGAGCCCGCTCGCCACCATCTAA
GCCACAGCCTTGCCAGGAGTCCACCACACCACCAAGTTAACCTGAGCCTAAACCGCAACAGCATACA
CCACAGCCTCTCCGCTCCAGAGAAACCTGCCTCATCTCCGTTCTGGTACCCATGGAACCTAAACCAA
TTCTTCCAAGCAGAAAGCAGCAGTGCCTCTGACCTATGTGGCCCCACGCTCCTGTCTGTGCAGCCTG
CCCGGGCAGCTCCGCTTGTGGCATCGTCTGGGCTTTGCCATAGCCGAATCTTTGATGTCTCCTACCT
CGGGACTGGTCATCCATGCCAGGGAGAGGAGTCCCAAACCTTCTCACCTTCTACAGGTTCCGCATGAACG
GGAGGGCTAAAAGCCTTGGGTG

ACGCGTACGCGGCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR201446 representing NM_027029
Red=Cloning site Green=Tags(s)

MKKVKKKSDSRRRRNSISPQTSSDSSQQPSSETPPSCPEPASPPSKPQPCQESTTPHQVNPEPKPQQHT
PQPLPPEKPASSPFLVPMEPKPILPSRKAAPLTYVAPRSCSACACPGSSACWHRLGLCHSRIFDVLLP
RDWSSMPGRGVPNLLTFYRFMNGRAKPPWV

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN: NM_027029

ORF Size: 510 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_027029.1</u> , <u>NM_027029.2</u> , <u>NP_081305.1</u>
RefSeq Size:	751 bp
RefSeq ORF:	513 bp
Locus ID:	70060
UniProt ID:	<u>Q9D9T6</u>
Cytogenetics:	1 C5
MW:	19.3 kDa