

Product datasheet for MR201455

Myl9 (NM_172118) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Myl9 (NM_172118) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: Myl9
Synonyms: AI327049; MLC20; Mylc2c; RLC-C
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >MR201455 representing NM_172118
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCCGCATCGCC

ATGTCGAGCAAGAGAGCCAAGGCCAAGACCACCAAGAAGAGGCCCCAGAGGGCTACGTCCAATGTCTTCG
 CCATGTTTGACCAGTCCCAGATCCAGGAGTTAAGGAGGCCTTCAACATGATTGATCAGAACCGAGATGG
 CTTTCATTGATAAGGAGGACCTGCACGACATGCTGGCCTCTCTGGGAAGAACCCACAGACGAGTATCTG
 GAGGGCATGATGAACGAGGCGCCAGGGCCATCAACTTACCATGTTCTCACAATGTTTGGGAGAAGC
 TGAACGGCACAGACCCCGAGGATGTGATCCGCAATGCCTTGCCTTGTGATGAGGAGGCCTCAGGCTT
 CATCCACGAGGACCACCTGAGGGAGCTGCTCACCACCATGGGCGACCGATTACCGGATGAGGAGGTGGAC
 GAGATGTACCGCGAGGCACCCATTGATAAGAAGGGCAACTTCAACTACGTGGAGTTCACCTCGCATCTCA
 AACACGGCGCCAAGGACAAGGACGAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR201455 representing NM_172118
 Red=Cloning site Green=Tags(s)

MSSKRAKAKTTKRPQRATSNVAMFDQSQIQEFKEAFNMIDQNRDGFIDKEDLHDMLASLGNPTDEYL
 EGMNEAPGPINFMTFLTMFGEKLNQDPEDVIRNAFACFDEEASGFIHEDHLRELLTTMGRDFTDEEVD
 EMYREAPIDKKGFNFYVEFTRILKHGAKDKDD

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Chromatograms: https://cdn.origene.com/chromatograms/mm9037_h05.zip



[View online »](#)

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_172118

ORF Size: 516 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_172118.1](#), [NP_742116.1](#)

RefSeq Size: 1083 bp

RefSeq ORF: 519 bp

Locus ID: 98932


UniProt ID: [Q9CQ19](#)

Cytogenetics: 2 H1

MW: 20.3 kDa

Gene Summary: Myosin regulatory subunit that plays an important role in regulation of both smooth muscle and nonmuscle cell contractile activity via its phosphorylation. Implicated in cytokinesis, receptor capping, and cell locomotion.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR201455