

Product datasheet for **MR201462L1V**

Cox4i2 (NM_053091) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Cox4i2 (NM_053091) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Cox4i2
Synonyms:	Cox4b; CoxIV-2
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_053091
ORF Size:	516 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR201462).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_053091.2 , NP_444321.1
RefSeq Size:	727 bp
RefSeq ORF:	519 bp
Locus ID:	84682
UniProt ID:	Q91W29
Cytogenetics:	2 H1
Gene Summary:	This protein is one of the nuclear-coded polypeptide chains of cytochrome c oxidase, the terminal oxidase in mitochondrial electron transport.[UniProtKB/Swiss-Prot Function]



[View online »](#)